




Nuvo Vu LED Surgical Light System



*Cross **F**ocus **T**echnology **S**ystem*




Cross Focus Technology enables the LED's to move independently of each other while adjusting the pattern size from 9.5" – 14" (229 mm – 356 mm). CFT allows the Vu to provide the same output on the large pattern as it does on the small pattern.



Ambient Light for MIS cases




Low profile for ease of positioning



Sterile handle control



<u>Intensity</u>	<u>160,000 lux</u>
<u>Color Temperature</u>	<u>4300°k</u>
<u>CRI</u>	<u>95</u>
<u>R9 Value</u>	<u>95</u>
<u>Pattern Size</u>	<u>9.5" – 14"</u>
	<u>229 – 356 mm</u>
<u>Depth of Illumination</u>	<u>32.3" (820.4 mm)</u>
<u>Intensity Settings</u>	<u>5</u>
<u>Ambient Light Mode</u>	<u>Yes</u>
<u>Controls on sterile handle</u>	<u>Yes</u>
<u>Life Expectancy of LED's</u>	<u>50,000 hours</u>
<u>Voltage</u>	<u>100 - 240 VAC 50/60 Hz</u>



Specifications subject to change without notice